



MXDPP-50 (Front)



MXDPP-50 (Back)

The MXDPP-50 digital pulse processor (DPP) is ideal for analytical x-ray and gamma-ray instruments such as; SDD, Si-PIN, CaTe, Si(Li), and Germanium detectors. The MXDPP-50 is designed to digitize detector output signals, achieving high throughput, with excellent pile-up rejection. The MXDPP-50 box includes the digital pulse processor, the detector power supply including high voltage bias, and the automatic temperature controller for Moxtek detectors.

### Features

### Benefits

|  |  |
|--|--|
| Analog-to-Digital Converter (ADC) system (50 MHz, 14 bit)                        | Fast processing speeds                               |
| Two independently configurable fast channels for pileup rejection                | Best pile-up rejection for wide range of energies    |
| Pulse / Continuous feedback  | Compatible with various detector types               |
| Positive / Negative step   |  |
| Configurable: HV bias, ramp polarity, temperature set point                      | Change detector settings without changing jumpers    |
| USB2 Interface   | High speed communication with Personal Computer (PC) |
| Certified Microsoft Windows Drivers XP, Vista, Windows 7 & 8 (32 bit and 64 bit) |  |


### Applications

- Portable XRF applications
- X-ray gamma ray detectors
- Process control
- Scientific research
- Nuclear monitoring

### Standard Package Includes

- MXDPP-50 Digital Pulse Processor
- Detector power supply with HV bias
- 5 VDC power supply brick
- XPIN detector cable (LEMO)
- XPIN signal cable (BNC)
- USB cable
- Software CD
  - Product manual (PDF)
  - X-spectrum DX software
  - Software development kit

### FP Software Sold Separately



Example MXDPP-50 Spectrum



## MXDPP-50 Product Specifications

### **Inputs:**

Detector signal input range -6 V to +6 V

- Software selectable signal polarity
- 2 Auxiliary software monitorable digital inputs

### **Multi Channel Analyzer(MCA):**

4G counts / channel maximum

Multi preset functions

- Real Time
- Live Time
- Total Counts
- Peak Counts

4096 channels

Live-time and Real-time acquisition (1 msec to 49 days)

### **Integrated non-volatile EEPROM Memory:**

For persistent parameter settings between power cycles

### **Outputs:**

Selectable analog output signal for monitoring filter outputs with oscilloscope or input into external MCA

2 Auxiliary digital outputs

8 Single Channel Analyzer (SCA) outputs with software selectable windows

- Selectable dead time corrected frequency output or pulse per event output

### **Pulse Processing:**

Trapezoidal shaping

Slow channel and two fast channels for pile-up rejection

Programmable peaking time (flat top) for all channels

- Slow 0.08 to 81 $\mu$ s
- Fast 0.08 to 19.8 $\mu$ s

Offset function for adjustment of 0 channel to 0eV

Adjustable holding time for all channels

- Slow 0.08 to 19.8  $\mu$ s
- Fast 0 to 19.8 $\mu$ s (typically 0)

Digital fine gain setting for adjustment of eV/channel with high resolution

### **Preprocessor:**

50MHz 14Bit ADC

Adjustable coarse gain for keeping same signal-to-noise ratio level at ADC

Adjustable differentiator time constant for optimal performance

### **Communication:**

Virtual Com Port (VCP) and Direct (D2XX) drivers available

USB 2.0 high speed communication

## Power Supply Specifications

### **Power Supply:**

+5V  $\pm$ 10% 1A Max (Includes power for detector)

Integrated detector power supply

- $\pm$ 5V for XT style detectors,  $\pm$ 9V for BT style detectors
- Software adjustable detector bias supply
  - -250Vdc to +250Vdc

### **Environmental:**

0 to 60°C Operating Temperature

EAR99 Export Classification

-40 to 80°C, 10 to 90% humidity non-condensing

Storage Temperature

RoHS Compliant

### **Integrated Detector Temperature Controller (TC):**

Software selectable TC mode

- Detector mode for use with detectors with integrated temperature controllers
- Box mode for use with detectors without integrated temperature controllers

Software selectable temperature setting

- 0°C to -100°C

Software monitor for Detector temperature, TEC voltage, Ready Signal, DPP temperature



## Free Evaluation Software (X-SpectrumDX) Specifications

### Read/Write DPP Parameters:

- Factory Default DPP configuration files for easy operation
- Fine Control of all DPP parameters for advanced users

### Read/Write Detector Settings:

- Predefined detector settings for Moxtek detectors
- Monitor and display detector temperature, TEC voltage

### Save Data:

- Save raw data to .csv file
- Save full data including ROIs, calibration, X-Ray lines to .csv file
- Save spectrum as image (.bmp, .gif, .jpg, .png, .tiff)

### Software Development Kit (SDK):

- Example Code
  - Visual Basic 2010
  - C#.NET 2010
  - LabView

### Intuitive GUI Interface:

- Linear/Logarithmic display
- Spectrum controls using mouse
- Integrated help for DPP parameters
- Copy & Paste spectrum data to/from excel using windows clipboard
- Customizable colors

### Basic XRF Analysis:

- Start/Stop acquisition controls
- Download spectra
- Automatic peak search
- 2 peak linear calibration
- Integrated X-Ray Line display
- Example spectrums

## Connector Specifications


### Connections:

Detector:


- Detector Power connector: Lemo, EPL.1S.306.HLN
- Detector Signal (BNC) connector: Amphenol, 31-5538-10RX

Auxiliary: Samtec, T2M-110-01-L-D-RA


Communications: USB type B connector



## DPP Signal Diagram and Spectrum




DPP Signal Example




DPP Spectrum Example



# Block Diagram of MXDPP-50 and XRF System



## MXDPP-50 Box Dimensions and Pinouts



| Detector Power Connector (Lemo*) |                   |                 |
|----------------------------------|-------------------|-----------------|
| Pin ID                           | Description       | Set             |
| 1                                | Temp / Ready      | Input           |
| 2                                | High Voltage Bias | -250 to +250Vdc |
| 3                                | Preamp Power -    | -9V Output      |
| 4                                | Preamp Power +    | +9V Output      |
| 5                                | Temp GND / TEC -  | TC GND / Return |
| 6                                | TC Power / TEC +  | Output          |
| Shield                           | Detector Ground   | GND             |

Table 1 Lemo Power Connector Pinout

| Signal Connector (BNC) |                 |
|------------------------|-----------------|
| Contact                | Description     |
| Pin                    | Detector Signal |
| Shield                 | Signal GND      |

Table 2 Signal Connector Pinout

Note: Dimensions are in (inches) and [millimeters]



452 West 1260 North / Orem, UT 84057  
 Phone 801.225.0930 / Fax 801.221.1121  
[www.moxtek.com](http://www.moxtek.com)  
[info@moxtek.com](mailto:info@moxtek.com)